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FILE OPTIONS
SAVE LOCATIONS
AUTO START
DEVICE CONFIGURATION
INTEGRATION

FILE PREFIX | `opensignals`

DEVICES | ☒

FILE FORMAT | ☐ .txt | ☒ .h5 | ☐ .edf

CONVERTED VALUES | ☐

The filename is generated with the FILE PREFIX followed by the name of each acquired device (if the DEVICES flag is enabled) and finished with the timestamp.

CONVERTED VALUES option will create a separate file containing the raw data converted to the corresponding physical units, as established by the sensor transfer function.

☐ CONTINUOUS MODE  
When this option is turned on, whenever the connection to a device is lost, the software will automatically try to reconnect and resume data acquisition

NOTICE

Some plugins require the .h5 in order to save their data. Be aware that if you choose to save using the .txt file format, that data will be lost.

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FILE OPTIONS
SAVE LOCATIONS
AUTO START
DEVICE CONFIGURATION
INTEGRATION

☐ AUTO-SAVE  
Turning on this option will allow you to save manually after each acquisition session.

SAVE LOCATION | `/Users/narimato@almare/Do` | 

REPOVIZZ | `[✓] repovizz` | ☐

DROPBOX |  | ☐

GOOGLE DRIVE |  | ☐

☐ CONTINUOUS MODE  
When this option is turned on, whenever the connection to a device is lost, the software will automatically try to reconnect and resume data acquisition

NOTICE

Some plugins require the .h5 in order to save their data. Be aware that if you choose to save using the .txt file format, that data will be lost.

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[1] BITalino (r)evolution PsychoBIT Bundle Data Sheet,

<https://www.creact.co.jp/measure/bio/bitalino/file/PsychoBIT-Bundle-Datasheet.pdf>.

[2] OpenSignals (r)evolution (Download),

<https://support.pluxbiosignals.com/knowledge-base/introducing-opensignals-revolution/>.